## RAPID SARS-COV-2 ANTIGEN TEST CARD

### INSTRUCTION GUIDE FOR ANTERIOR NASAL SWAB SPECIMENS

#### For self-testing

0123

REF	1N40C5-2	For 1 Test/Box
REF	1N40C5-4	For 5 Tests/Box
REF	1N40C5-6	For 20 Tests/Box

Please follow the instruction leaflet carefully.

**INTENDED USE**  
Rapid SARS-CoV-2 Antigen Test Card is an immunochromatography based one step in vitro test. It is designed for the rapid qualitative determination of SARS-CoV-2 virus antigen in anterior nasal swabs from individuals with or without symptoms or other epidemiological reasons to suspect COVID-19. Rapid SARS-CoV-2 Antigen Test Card does not detect the SARS-CoV-2 nucleocapsid protein (N protein). Rapid SARS-CoV-2 Antigen Test Card shall not be used as sole basis to diagnose or exclude SARS-CoV-2 infection. Children under 14 years of age should be assisted by an adult.

**SUMMARY**  
The novel coronavirus belongs to the *Beta* genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection, asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

Components	For 1 Test/Box	For 5 Tests/Box	For 20 Tests/Box
Rapid SARS-CoV-2 Antigen Test Card (sealed foil pouch)	1	5	20
Sterile swab	1	5	20
Extraction tube	1	5	20
Sample extraction buffer	1	5	20
Instructions for use (this leaflet)	1	1	1
Tube stand	1 (packaging)	1	1

**PERFORMANCE (SENSITIVITY AND SPECIFICITY)**  
Rapid SARS-CoV-2 Antigen Test Card was compared to the confirmed clinical diagnosis. The study involved 1063 nasal samples. The test results are summarized below:

Evaluated Reagent Results	RT-PCR Results		Total
	Positive (+)	Negative (-)	
Positive (+)	425	1	426
Negative (-)	10	627	637
Total	435	628	1063

Sensitivity (PPA) = 425/435\*100% = 97.70% (95% CI: 96.29% - 99.11%)  
Specificity (NPA) = 627/628\*100% = 99.84% (95% CI: 99.53% - 99.99%)  
Accuracy (OPA) = 1052/1063\*100% = 98.97% (95% CI: 98.36% - 99.57%)  
A feasibility study demonstrated that:  
- 99.84% of non-professionals carried out the test without requiring assistance  
- 99.82% of the different types of results were interpreted correctly

**INTERFERENCES**  
None of the following substances at the tested concentration showed any interference with the test:

Whole Blood:	Alkalot: 10%	Mucin: 2%	Oxymetazolin: 15%
Phenylephrine: 15%	Tobramycin: 0,0004%	Oxymetazolin: 15%	Oxymetazolin: 15%
Menthol: 0,15%	Cromolyn: 0,15%	Benzocaine: 0,15%	Benzocaine: 0,15%
Fluticasone Propionate: 5%	Mupirocin: 0,25%	Zincam Nasal Spray: 5%	Zincam Nasal Spray: 5%
Oseltamivir Phosphate: 0,5%	Sodium chloride: 5%	Human Anti-mouse Antibody (HAMA): 60 ng/mL	Human Anti-mouse Antibody (HAMA): 60 ng/mL

**CROSS REACTIVITY**  
Cross-reactivity of the Test Device was evaluated by testing viruses and other microorganisms. The final test concentrations of viruses and other microorganisms are documented in the Cross-Reactvity Study. The following viruses and other microorganisms except the Human SARS-coronavirus have no effect on the test results of the Test Device: Human coronavirus 229E, Human coronavirus OC43, Human coronavirus NL63, Human coronavirus HKU1, MERS coronavirus, Parainfluenza Virus 1-4, Enterovirus EV71, Respiratory syncytial virus, Rhinovirus, Influenza A virus (H1N1 and H3N2), Influenza B virus (Yamagata and Victoria), Adenovirus (HUMAN METAPNEUMOVIRUS (hMPV)), Staphylococcus epidermidis, Chlamydia pneumoniae, Streptococcus pneumoniae, Streptococcus pyogenes, Bordetella pertussis, Mycobacterium tuberculosis, Legionella pneumophila, Mycoplasma pneumoniae, Haemophilus influenzae, Candida albicans, Staphylococcus aureus, Pseudomonas aeruginosa, Escherichia coli, Pneumocystis jirovecii (PJP) and pooled human nasal secretions.

**IMPORTANT INFORMATION BEFORE THE EXECUTION**  
1. Read this instruction guide carefully.  
2. Do not use the product beyond the expiration date.  
3. Do not use the product if the pouch is damaged or the seal is broken.  
4. Store the test device at 4 to 30 °C in the original sealed pouch. Do Not Freeze.  
5. The product should be stored at room temperature (15° to 30 °C). If the product has been stored in a cool area (less than 15 °C), leave it at normal room temperature for 30 minutes before using.

6. Handle all specimens as potentially infectious.  
7. Inadequate or inappropriate specimen collection, storage, and transport may yield inaccurate test results.  
8. Use the swabs included in the test kit to ensure optimal performance of the test.  
9. Correct specimen collection is the most important step in the procedure. Make sure to collect enough specimen material (nasal secretion) with the swab, especially for anterior nasal sampling.  
10. Blow the nose several times before collecting specimen.  
11. The specimens should be tested as soon as possible after collection.  
12. Apply the drops of test specimen only to the specimen well (S).  
13. Too many or too few drops of extraction solution can lead to an invalid or incorrect test result.  
14. When used as intended, there should be no any contact with the extraction buffer. In case of contact with skin, eyes, mouth or other parts, rinse with clear water. If an irritation persists, consult a medical professional.

15. Children under 14 years of age should be assisted by an adult.

**LIMITATIONS**  
1. The test is to be used exclusively for the qualitative detection of SARS-CoV-2 viral antigen in anterior nasal swab specimens. The exact concentration of SARS-CoV-2 viral antigen in the Cross-Reactvity Study. The following viruses and other microorganisms except the Human SARS-coronavirus have no effect on the test results of the Test Device: Human coronavirus 229E, Human coronavirus OC43, Human coronavirus NL63, Human coronavirus HKU1, MERS coronavirus, Parainfluenza Virus 1-4, Enterovirus EV71, Respiratory syncytial virus, Rhinovirus, Influenza A virus (H1N1 and H3N2), Influenza B virus (Yamagata and Victoria), Adenovirus (HUMAN METAPNEUMOVIRUS (hMPV)), Staphylococcus epidermidis, Chlamydia pneumoniae, Streptococcus pneumoniae, Streptococcus pyogenes, Bordetella pertussis, Mycobacterium tuberculosis, Legionella pneumophila, Mycoplasma pneumoniae, Haemophilus influenzae, Candida albicans, Staphylococcus aureus, Pseudomonas aeruginosa, Escherichia coli, Pneumocystis jirovecii pneumonia (PJP) and pooled human nasal secretions.

2. The SARS-CoV-2 rapid antigen test can identify both viable and non-viable SARS-CoV-2 material. The performance of the SARS-CoV-2 rapid test is dependent on viral load and may not correlate with other diagnostic methods performed on the same specimen.  
3. Users should test specimens as soon as possible after specimen collection and within two hours of specimen collection.  
4. Sensitivity for nasal or oropharyngeal swabs may be lower than nasopharyngeal swabs. It is recommended to use the nasopharyngeal swab specimens by healthcare professionals.  
5. Monoclonal antibodies may fail to detect, or detect with less sensitivity, SARS-CoV-2 viruses that have undergone minor amino acid changes in the target epitope region.  
6. The amount of antigen in a sample may decrease as the duration of illness increases. Specimens collected after day 5-7 of illness are more likely to be tested negative compared to a RT-PCR assay.

7. The kit was validated with the associated swabs. Use of alternative swabs for identification/confirmation of tissue culture isolates and should not be used in this capacity.  
8. The validity of Rapid SARS-CoV-2 Antigen Test Card has not been proven for identification/confirmation of tissue culture isolates and should not be used in this capacity.  
9. Positive test results do not rule out co-infections with other pathogens. Positive results may occur in cases of infection with SARS-CoV-2.  
10. Clear, clean and dry a flat surface.  
11. Check the test kit contents. Make sure that nothing is damaged or broken.  
12. Tear off the seal.  
13. Blow your nose several times before collecting specimen.  
14. Wash hands.

**DISPOSAL**  
The test kit may be disposed of with normal household waste in accordance with the applicable local regulations.

**PROCEDURE**  
This test is suitable for people of all ages. The recommended operators are aging from 14 to 90. Children under 14 years of age should be tested by an adult. Do not continue the test if the child feels any pain.

1. Rotate the lid of sample extraction buffer bottle. **Caution: Open it away from your face and be careful not to spill any of the liquid.**

2. Squeeze the air extraction buffer out of the bottle into the extraction tube. **Caution: Avoid touching the bottle against the tube.**

3. Find the swab in the sealed wrapper in front of you. Identify the soft, fabric tip of the swab.

4. Peel open the swab packaging and gently take out the swab. **Caution: Never touch the soft, fabric tip of the swab with your hands.**

5. Carefully insert swab into one nostril. The swab tip should be inserted no less than 2.5 cm (1 inch) from the edge of the nostril. Roll swab 3-4 times along the mucosa inside the nostril. Leave swab in the nostril for several seconds. Using the same swab, repeat this process for the other nostril. Withdraw swab from the nasal cavity. **Caution: This may feel uncomfortable. Do not insert the swab any deeper if you feel strong resistance or pain.**

6. Place swab into extraction tube. Roll swab three to five (3-5) times. **Leave swab in extraction buffer for 1 minute.**

7. Pinch extraction tube with fingers and remove the solution from swab as much as possible.

8. Bring the kit components to room temperature before testing. Open the pouch and remove the card. Place the card on a flat and level surface. **Caution: Once opened, the test card must be used immediately.**

9. Invert the extraction tube and add 3 drops (about 75 µL) of test specimen into the specimen well (S) by gently squeezing the extraction tube. **Caution: The formation of air bubbles in the specimen well (S) must be avoided.**

10. Read the results at 15-20 minutes. **Caution: Results after 20 minutes may not be accurate. The used device may be disposed of with normal household waste in accordance with the applicable local regulations.**

11. Interpretation of results: **Positive:** If two colored bands appear with one colored band in the Control Zone (C) and another in the Test Zone (T) within 15-20 minutes, the test result is positive. **Caution: No matter how faint the colored band is in the Test Zone (T), the result should be considered as positive.**

**Negative:** If one colored band appears in the Control Zone (C) and no colored band appears in the Test Zone (T) within 15-20 minutes, the test result is negative.

**Invalid:** If no color line appears in the control area (C) within 15-20 minutes, the test is invalid. Repeat the test with a new test card.

**QUALITY CONTROL**  
The control line is an integrated reagent and is used to control the procedure. The control line appears when the test has been performed correctly and the reagents are reactive.

**FREQUENTLY ASKED QUESTIONS (FAQ)**  
1. How does the detection work?  
The reaction of the SARS-CoV-2 virus reacts with the stripe-like coating of the test line and, if present, results in a color change, i.e. a red line appears. Therefore, if the sample does not contain any viral proteins or antigens, there will be no red test line (T).

2. When should I can I test myself?  
You can test yourself whether you have symptoms or not. Studies show that earlier testing within the first 4 days of illness typically means a higher viral load, which is easier to detect. Since the test result is a snapshot valid for that point in time, testing should be repeated as recommended by local authorities.

3. What can affect my test result? What should I pay attention to?  
Be sure to follow your nose multiple times before collecting the specimen.  
Be sure to visibly collect sample material (nasal secretions).  
Perform the test immediately after taking the sample.  
Follow the instructions for use carefully.  
Apply the drops of extraction solution only to the sample well (S).  
Too many or too few drops of extraction solution can lead to an invalid or incorrect test result.

4. The test strip is clearly discolored or smudged? What is the reason for this?  
Please note that the test card should not be used with more than 3 drops of sample, as the liquid absorption of the test strip is naturally limited. If the control line does not appear or the test strip is badly smudged or discolored, making it unreadable, please repeat the test according to the instructions.  
I have taken the test, but I don't see a control line (C). What should I do?  
Your test result is invalid. Observe the answer to question 4 and repeat the test according to the instructions for use.

5. I am unsure about reading the result. What should I do?  
For the result to be positive, 2 straight horizontal lines must be clearly visible with the full width of the cassette. If you are still unsure about the results, contact the nearest health facility according to the recommendations of your local authorities.

6. My result is positive. What should I do?  
Your result is positive and the test kit clearly indicates the control line as well as the test line, you should contact the nearest medical facility as recommended by your local authorities. Your result may be double-checked at any authority or facility, which will explain the appropriate next steps.

7. My result is negative. What should I do?  
If the test kit only clearly shows the control line, this may mean that you are negative or that the viral load is too low to be detected. If you experience symptoms (fever, cough, mild loss of sense of smell or taste, etc.), test and consult your primary care physician, or the nearest health care facility as recommended by your local authorities.  
If you are not sure, you can repeat the test.

8. How can I dispose of the product?  
The test kit may be disposed of with normal household waste in accordance with the applicable local regulations.

9. How can I dispose of the product?  
The test kit may be disposed of with normal household waste in accordance with the applicable local regulations.

Accessories	Manufacturer	EC-Representative	CE-Mark
Swab A	Jiangsu Changfeng Medical Industry Co., Ltd Touqiao Town, Guangling District Yangzhou 225109 Jiangsu P.R. China	Lins Service & Consulting GmbH Obere Seegasse 34/2,69124 Heidelberg Germany	117 acc. 93/42/EEC
Swab B	Goodwood Medical Care Ltd. 1-2 Floor, 3-919 Yongzheng Street, Jinzhou District Dalian 116100 Liaoning China	CMC Medical Devices & Drugs S.L. C/ Horacio Lengua No18, CP 29006, Málaga, Spain	197 acc. 93/42/EEC
Swab C	Zhejiang Gongdong Medical Technology Co., Ltd. No. 10 Beiyuan Ave., Huangyan 318020 Taizhou, Zhejiang, P.R.China	Shanghai International Holding Corp. GmbH (Europe) Eiffelstrasse 80, 20537 Hamburg, Germany	123 acc. 93/42/EEC
Swab D	Jiangsu Hanheng Medical Technology Co., Ltd. 16-84, #1 North Qinyang Road, Tiaanning District, 213017, Changzhou, Jiangsu, China	Luxus Lebenswelt GmbH Kochstr. 1, 47877, Willich, Germany	197 acc. 93/42/EEC

IVD	In Vitro Diagnostics Use	See Instructions for Use	Expiry Date
Tests per Kit	Keep dry	LOT	Batch Number
Authorized Representative	Keep away from sunlight	Manufacturer	Manufacturer
Do not reuse	Do not use if package is damaged	Store between 4-30 °C	Warning, please refer to the instruction

CE Mark	REF	Catalogue Number	Warning, please refer to the instruction
0123	1N40C5-2	1N40C5-2	1N40C5-2
1N40C5-4	1N40C5-4	1N40C5-4	1N40C5-4
1N40C5-6	1N40C5-6	1N40C5-6	1N40C5-6

STERILE/EO	Sterilized using ethylene oxide	Authorized Representative:	Lotus N. B.V. Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.
Manufacturer:	Xiamen Boson Biotech Co., Ltd. 90-94 Tianfeng Road, Jimei North Industrial Park, Xiamen, Fujian, 361021, P.R.China.	Authorized Representative:	Lotus N. B.V. Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.
Version 6.0	Date: August 15, 2021	Version 6.0	Date: August 15, 2021

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